

CONDUCTOMETRY -INTRODUCTION, CONDUCTIVITY CELL

BY, LIPSA SAMAL ASST. PROF (PA & QA), SPLS, CUTM



INTODUCTION:

- It is an electrochemical method of analysis concerned with electrical conductance through an electrolyte solution . (or) It is defined has as determination or measurement of the electrical conductance of an electrolyte solution by means of a conductometer.
- electric conductivity of an electrolyte solution depends on :
 - **1.** Type of ions (cations, anions, singly or doubly charged
 - 2. Concentration of ions
 - 3. Temperature
 - 4. Mobility of ions

PRINCIPLE:

enturio

Based on the conductance of electrical current through electrolyte Empowering Communities... solutions similar to metallic conductors.

> The electric conductance in accordance with ohms law which states that the stength of current(i) passing through conductor is **directly** porportional to potential difference & inversely to resistance. i = V/R.

Resistance (r):

Is a measure of the conductors opposition to the flow of electric charge R =1/G

<u>Specific resistance:(p)</u>

Is resistance offered by a conductor of unit length and having unit **cross section** $\mathbf{R} \propto \mathbf{l}/\mathbf{A}$



Conductance:(G)

Conturion UNIVERSITY with which current flows per unit area of conductor per unit Shaping Lives tential applied & is reciprocal to resistance(R)

Specific conductance (K):

conductance of the body of uniform length(l)& uniform area cross section

Molar conductance:

Conductance of a solution containing 1 mole of the solute in 1000 cm of the solution which placed between two parallel electrodes which are 1 cm apart

Equivalent conductance:

specific conductance of the solution containing 1gm equivalent of solute in 1000cm3 of solution.



CONDUCTIVITY CELLS:-

- Made of pyrex or quartz and are fitted with two platinum electrodes.
- Should be placed in vessel containing water to maintain constant temperature
- Types : 1. Wide mouthed cell 2. Cell for reactions producing precipitation
- <u>Wide mouthed cell:-</u>
- Measurement of low conductance
- Wide mouthed fitted with an ebonite cover which has provisions for platinum electrodes and burettes



Cell for reactions producing ppts:

- Mainly used for ppt reactions
- Also wide mouthed fitted with ebonite cover which has provisions for burette ,electrode , stirrer
- Stirrer may be mechanical or magnetic



ELECTRODES:

- Platinum sheets, each of 1 cm2 are fixed at distance of 1 cm
- The surface is coated with platinum black to avoid polarization effects and increase effective surface area.
- Platinisation of electrodes is done by coating solution of 3% chlorplatinic acid and lead acetate on it to get uniform coating
- **Electrodes usage depends on conductivity and concentration**
- If conc is low then electrodes should be largely and closely packed



A DIAGRAM OF A TYPICAL ELECTROCHEMICAL CELL





DIAGRAM OF A TYPICAL ELECTRODE



THANK YOU