Experiment – 11

<u>Aim:-</u> Preparation of glycerogelatin base suppositories

Rx

Gelatin 14 g Glycerin 70 g Water QS to 100g

Equipments Required:-

Mortar and pestle, Porcelain crucible, water-bath, 2 g x 6 mold

Chemicals:

Gelatin 84g, Glycerin 420 g, water 96 ml.

Procedure:

- 1. Calculate for 6 supp., that is the amount of Gelatin = $6 \times 14 = 84 \text{ g}$.
- 2. Weigh the required amount of gelatin and soak it in enough water until thoroughly softened.
- 3. Put the soaked gelatin on the calculated amount of glycerin previously weighed in a tared crucible and on a water bath until gelatin dissolves and a constant weight is obtained.
- 4. During evaporation the liquid mass should be only gently stirred, rapid stirring produce air bubbles which may appear in the finished suppositories.
- 5. Remove any skin formed on the surface before pouring.
- 6. Pour the mass while still hot, into the holes lubricated with liquid paraffin, do not let the melted mass to overflow.