



**Centurion**  
**UNIVERSITY**  
*Shaping Lives...  
Empowering Communities...*

# **PREPARATION AND STANDARDIZATION OF CERIC AMMONIUM SULPHATE**

**BY, LIPSA SAMAL**  
**ASST. PROF (PA & QA), SPLS, CUTM**



**Centurion**  
**UNIVERSITY**

*Shaping Lives...  
Empowering Communities...*

## **STANDARDIZATION OF CERIC AMMONIUM SULPHATE**

**Ceric ammonium sulphate solution is titrated with arsenic trioxide (primary standard) in presence of sulfuric and osmic acid utilizing the ferroin sulphate as indicator as the same method as oxallic acid.. End point is visually detected by change in color from pink to very pale blue.**



**Centurion**  
**UNIVERSITY**

*Shaping Lives...*  
*Empowering Communities...*

**THANK YOU**